

ONE MINUTE MEDIC

EL X RAPID

PACK Works with just 60 seconds compression¹

FIGHT Remain engaged - suppress the threat

LIVE Stops life-threatening bleeding fast²

Contact us at info@celoxonline.com

CELOX RAPID STOP BLEEDING. SAVE LIVES.

THE NEW STANDARD FOR HEMORRHAGE CONTROL

Time is critical when treating severe injuries in the field and saving time increases the chances of survival for the casualty in an emergency or hostile situation. To be effective, hemostats need to work fast.

The current generation of hemostatic agents, Quikclot Combat Gauze LE and standard Celox™ Gauze require at least 3 minutes of compression to work. Celox™ Rapid Gauze removes the delay by speeding up packing time and reducing compression time to just 60 seconds.



CUT TREATMENT TIME IN HALF

Because Celox[™] Rapid sticks to wet tissue, it only needs 60 seconds of compression with no loss of performance; What is more, Celox[™] Rapid also saves time packing the wound⁴, so the real saving is more than two minutes.

Ordering Information

Celox Rapid 3" x 5' Z-fold Gauze: FG08839011 NSN: 6510-01-623-9903 (single dressing) 6510-01-623-9456 (case of 100 dressings) Celox Rapid Ribbon 1" x 5' Z-fold gauze: FG08839041 NSN: 6510-01-623-9442 (single dressing) 6510-01-623-9730 (case of 150 dressings)

Biologistex, Inc.

171 Cooper Avenue, Suite 106, Tonawanda, NY 14150 Toll Free 866-992-3569

HOW CELOX™ RAPID WORKS

Celox[™] Rapid contains Chito-R[™] activated chitosan granules that stick to wet tissue in the wound^³.

- Celox[™] Rapid absorbs blood fluid and forms an adhesive gel that seals the wound to stop the flow of blood, speeding up compression time;
- Works on hypothermic blood;
- No minerals, chitosan is made up of glucosamine and N-acetyl glucosamine.

CELOX PROVEN PERFORMANCE

Celox[™] Gauze testing by Office of Naval Research on behalf of DoD⁵.



Celox™ Rapid uses the same chitosan technology as Celox Gauze, which obtained highest observed survival in US DoD testing.





Combat Gauze XL

Chitogauze™ Combat Gauze™

WHAT THIS MEANS

- Celox[™] Rapid has been independently shown to reduce blood loss by 60% compared to QuikClot Combat Gauze LE²
- Celox[™] Rapid stays in place during evacuation³
- Also available; Celox Rapid Ribbon, a compact carry solution for small narrow entry wounds.

Contact us at info@celoxonline.com

1. Product instructions for use. 2. Chitosan based advanced hemostatic dressing is associated with decreased blood loss in swine uncontrolled hemorrhage model. Kunio et al. Am J Surg 2013. 3. Mechanism of action of rapid-acting gauze hemostat. A Hoggarth et al. Poster at ATACCC 1, FL. 4. Reduced application time with a rapid-packing gauze hemostat. A Hoggarth et al. Poster at ATACCC 1, FL. 5. Rall JM et al. Comparison of novel hemostatic dressings with QuikClot combat gauze in a standardized swine model of uncontrolled hemorrhage. J Trauma 2013. Trauma Gauze result (50%) not shown in graphic - product discontinued. TM Trademark of Medtrade Products Ltd. *Trademark of Z-Medica Corporation. Image UK MOD Crown Copyright 2015